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membrane as MP changes. MP depolarization results in an increase in fluorescence intensity of cells at Ex/Em = 590/616nm. Conversely, hyperpolarization is indicated by a decrease in fluorescence. PMA-differentiated THP-1 cells were cultured in 35mm imaging dish with a concentration of 4×10^5 cells/ml. Cells were stained in 200 nM DiBAC4(5) in 5mM K⁺ buffer at 37°C for 30 min. Then the staining buffer was removed without washing and crystals (silica, Al(OH)₃, CPPD) or inhibitors (glyburide or CFTRinh-172) in fresh 5mM K⁺ buffer containing 200nM DiBAC4(5) was added to the cells. For high K⁺ condition, staining buffer was replaced with 100mM K⁺ buffer containing 200nM DiBAC4(5). We captured 5 different fields at each time-point on fluorescent microscope, and analyzed the fluorescence intensity of all the cells by ImageJ.

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
Fig. Depolarization or Hyperpolarization of Membrane Potential May Affect the Activation of NLRP3 Inflammasome. THP-1 cells were stained with DiBAC4(5) to check the effect of CFTRinh-172 on MP. Representative images and relative fluorescence unit (RFU) plotting are shown. Each point was mean $\pm$ SEM of cells from at least 12 fields in 2 independent experiments. Scale bar, 20 μ m.

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