**FM 1-43 (N-(3-Triethylammoniumpropyl)-4
-(4-(dibutylamino)styryl) pyridinium dibromide)
神经末梢荧光探针**

产品编号	产品名称	包装规格
NBS3210-1mg	FM 1-43 (N-(3-Triethylammoniumpropyl)-4 -(4-(dibutylamino)styryl) pyridinium dibromide)	1mg

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◇ FM 1-43 标记液: Add 1 μ l FM 1-43 stock solution to 1 mL HL-3 plus 90 mM KCl solution with calcium (final concentration of FM 1-43 is 4 μ M).

◇ Incubate the dissected larva in FM 1-43 labeling solution by replacing the HL-3 solution without calcium. Do not add the solution on top of the larva but gently pipet the solution on the side of the larva. Start a timer; once you add the labeling solution, exocytosis and endocytosis are induced, and FM 1-43 labeling ensues. The incubation time will depend on your experiment, but a 1-min incubation leads to robust labeling in wild-type animals). Since some mutants are temperature sensitive, it may be necessary to perform experiments at high temperature (see Note 2 if this is the case).

◇ Following stimulation, remove FM 1-43 labeling solution and wash five times over a total of 5-10 min with a generous amount of HL-3 solution without calcium to stop stimulating and to remove dye not internalized. Do not add the solution on top of the larva but gently pipet the solution on the side of the larva to avoid muscle damage. Gently perfuse the wash solution by pipeting up

and down.

◇Image labeled vesicles on a confocal microscope with an $\times 40$ water immersion lens and quantify intensity of labeling as described in Subheading 3.5.


Fig. FM 1–43 labeling of a third-instar larval fillet. (A)–(E) Third instar larval fillet incubated for 1 min in HL-3 with 4 μ M FM 1–43, 90 mM KCl, and 1.5 mM calcium, washed, and the FM 1–43 imaged (D). The same preparation was then briefly fixed with 3.7% formaldehyde and labeled with mouse anti-DLG antibodies (28; 4F3 monoclonal antibody; Developmental Studies Hybridoma bank) used at 1:50 (B), (C), and (E) to reveal synaptic boutons at the NMJs and Alexa 635-conjugated phalloidin (Invitrogen cat. no. A34054) used at 0.001 unit/ μ L (A) and (C) to reveal muscles. (A)–(C) two hemisegments labeled with phalloidin (A) and DLG (B), merged in (C). The dashed line in (A) shows the ventral midline, and muscles commonly used in the genetic analysis of synaptic function are marked. (D) FM 1–43 labeling (before fixation of the preparation) of the NMJ in the boxed area on muscle 4 shown in (C). (E) The same synapse as shown in (D), labeled with DLG (see B) following fixation of the preparation.

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